

<b>BRUKER AURORA M90, 810-MS and 820-MS</b>		
<b>Meinhard P/N°</b>	<b>OEM P/No.</b>	<b>Descriptions of parts</b>
TRP-50-A0.5		Meinhard Type A+ Concentric Glass Nebulizer, with PFA argon and sample connectors
TQP-50-A0.5		Meinhard Type A+ Concentric Quartz Nebulizer, with PFA argon and sample connectors
SiIQ+		High-Purity Quartz concentric nebulizer, with PFA fittings as standard
TR-50-K2	2010069500	K-Type Concentric Glass Nebulizer, 2 ml/min.
TR-30-K1	2010081600	K-High Flow, AX Conc. Glass Nebulizer, 1 ml/min.
TR-30-C2	2010076500	C-Type standard Concentric Glass Nebulizer, 2 ml/min.
PMM 4000		PEEK MiraMist Parallel-Path Nebulizer
PRL-01	8100-2051	'PEARL' Cleaning Device for concentric glass nebulizers - suits all types, all manufacturers
PRG-01		'PERGO' argon humidifier - keeps nebulizers running smoothly
PRG-01-MAPlus		'PERGO' argon high solids accessory including Meinhard Type A Plus high performance nebulizer
ML180031	0102-932	Quartz Double Pass Spray Chamber, ICPMS
ML180033	0102-933	Bruker Quartz Make-up IN-Port 2
ML180032	0102-1064	Quartz Universal Joint, Extension
ML18005S		Quartz Torch, Radial with Ball Joint
ML180030	2010100700	Quartz Torch, ICPMS, with Ball Joint, 2 quartz locating beads
ML180006	2010078300	Quartz Torch with 1.5mm-ID Injector
ML180011	2010079390	Bruker Quartz Torch, Demountable ICP-MS
ES-4398-1076		Peri-pump tubes for Sample, 2-tags, black-black, ID = 0,76 mm, pkt 12
ES-4398-1102		Peri-pump tubes for Sample, 2-tags, white-white, ID = 1,02 mm, pkt 12
ES-4398-1165		Peri-pump tubes for Drain, 2-tags, blue-blue, ID = 1,65 mm, pkt 12
ES-4397-1038		PVC 2-Stop Tubing, flared ends, ID=0.38 mm, Orange/Green, pkt 12
ES-4397-1051		PVC 2-Stop Tubing, flared ends, 0.51 mm, Orange/Yellow, pkt 12
Tube-PVCSD-51		PVC 2-Stop Tubing 0.051" Gray/Gray, pkt 12
Tube-PVCSD-65		PVC 2-Stop Tubing 0.065" Blue/Blue, pkt 12