

Components for Finnigan Element 1, 2 , XR and Neptune		
Meinhard #	OEM P/N*	Description
TQ-30-A1	1110720	Type A nebulizer, Quartz, 1 mL/min and 1 L/min @ 30 psi
TQP-50-A0.5		Type A-plus nebulizer, Quartz, low uptake, 1 L/min @ 50 psi
SilQ+		High-Purity Quartz concentric nebulizer, with PFA fittings as standard
TR-30-A1		Type A nebulizer, borosilicate glass, 1 mL/min and 1 L/min @ 30 psi
TR-50-A1		Type A nebulizer, borosilicate glass, 1 mL/min and 1 L/min @ 50 psi
TR-50-A0.2	1161870	Type A nebulizer, borosilicate glass, 0.2 mL/min and 1 L/min @ 50 psi
F0		PFA Capillary with encapsulated Micro Tubing. L=70cm, ID=0.25mm. For low uptakes.
F1		Detachable PFA capillary, for general application, L=70cm, ID=0.8mm.
F2		Standard PFA Integrated Capillary, for self-aspiration or pumped samples, L=70cm, ID=0,50mm.
FG		PFA Gas Quick Connect Fitting for Meinhard nebulizers. Zero contamination, easy to use.
PRL-01		'PEARL' Cleaning Device for concentric glass & quartz nebulizers - suits ALL types
PRG-01		'PERGO' argon humidifier, NO pressure vessel - based on membrane permeation
PRG-01-MAPlus		'PERGO' argon high solids accessory including Meinhard Type A Plus high performance nebulizer
ML185018	1120900	Quartz Spray Chamber, Element, Jacketed for Coolant, Threaded Port Connections <i>Cyclonic and mini-cyclonic quartz spray chambers are also available – please ask for details</i>
ML185026	ES-1001-0001	Torch Bonnet, quartz, L = 26 mm.
ML185025	1136500	Quartz Collar/Sleeve 60 mm, ES-1001-0002
ML185006		Quartz Torch, Element, 1.0 mm Injector
ML185007	ES-1002-0001	Quartz Torch, Element, 1.5 mm Injector
ML185008		Quartz Torch, Element, 2.0 mm Injector
ML185010	1091250	Quartz Torch, Element, Demountable Tube
ML185040	1130680	Quartz Injector with O-ring Ball Joint, 1.00 mm (does not include adapter), Thermo # 1161910
ML185041		Quartz Injector with O-ring Ball Joint, 1.50 mm (does not include adapter)
ML185042	1122800	Quartz Injector with O-ring Ball Joint, 1.80 mm (does not include adapter)
ML185043		Quartz Injector with O-ring Ball Joint, 2.00 mm (does not include adapter)
ML185033		Quartz Extension
ME10021000	ES-1002-1000	O-ring free, Precision Quartz Torch, Element, ESI
ME1002	ES-1002	Quartz Torch tube, Type CD
ES-1001-0004	1126640	Platinum torch guard shield, Length= 84 mm., Pt ring diam 20 mm
ES-4397-1038		PVC 2-Stop Tubing, flared ends, ID=0.38 mm, Orange/Green
ES-4397-1051		PVC 2-Stop Tubing, flared ends, 0.51 mm, Orange/Yellow
Tube-PVCSD-51		PVC 2-Stop Tubing 0.051" Gray/Gray
Tube-PVCSD-65		PVC 2-Stop Tubing 0.065" Blue/Blue