

GV-Micromass IsoProbe and Platform ICP-MS		
Meinhard P/N°	OEM P/No.	Part Descriptions
TQ-30-A2		Quartz concentric nebulizer, 2 mL/min, 1 L/min, 30 psi
TRP-50-A0.5		Meinhard Low-uptake Type-A+ Concentric Glass Nebulizer
TQP-50-A0.5		Meinhard Low-uptake Type-A+ Quartz Concentric Nebulizer
SilQ+		High-Purity Quartz concentric nebulizer, with PFA fittings as standard
TR-30-A1		Type A Concentric Glass Nebulizer, high sensitivity
TQ-30-A2		Quartz Type-A concentric nebulizer, 2 mL/min, 1 L/min, 30 psi
F0		PFA Capillary with encapsulated Micro Tubing. L=70cm, ID=0.25mm. For low uptakes.
F1		Detachable PFA capillary, for general application, L=70cm, ID=0.8mm.
F2		Standard PFA Integrated Capillary, for self-aspiration or pumped samples, L=70cm, ID=0,50mm.
FG		PFA Gas Quick Connect Fitting for Meinhard nebulizers. Zero contamination, easy to use.
PRL-01	8100-2051	'PEARL' Cleaning Device for concentric glass nebulizers - suits all types, all manufacturers
PRG-01		'PERGO' argon humidifier - keeps nebulizers running smoothly
PRG-01-MAPlus		'PERGO' argon high solids accessory including Meinhard Type A Plus high performance nebulizer
ML165016	6368303	Pyrex Spray Chamber, Double Pass
ML165015	6368308	Quartz Spray Chamber, Double Pass
ML165020	M695210BD1	Quartz Elbow, Ball/Socket
MM151030		Teflon Adapter (end cap), 6mm Nebulizer to 35mm Spray Chamber
ML165005	6368301	Quartz Torch, 1.5mm Injector
ML165008	6366307	Quartz Torch, Demountable for ISIS
ML165025	M695531A	Quartz Bonnet
EP-1202-1802		Semi-Demountable Torch, O-ring free, with Sapphire Injector and spray chamber adapter
EP-1203-0182		Sapphire injector for O-ring free Semi-Demountable Torch, with spray chamber adapter
EP-1202-0002		O-ring free Semi-Demountable Torch with PFA base, without injector
Tube-PVCSD-40-3		PVC 3-Stop Tubing 0.040" (1.02mm) White/White
Tube-PVCSD-075-3		PVC 3-Stop Tubing 0.0075" (0.19 mm) Orange/Red
Tube-PVCSD-075-3F		PVC 3-Stop Tubing 0.0075" (0.19 mm) Orange/Red, Flared both ends
Tube-PVC-10-FA		PVC 2-stop (72 mm) Tubing, 0.010" (0.254 mm), Orange/Blue, Flared Ends, 19 cm
Tube-PVC-075-FA		PVC 2-stop (72 mm) Tubing, 0.0075" (0.19 mm), Orange/Red, Flared Ends, 19 cm