

Meinhard parts for Horiba - Jobin-Yvon ICP instruments		
Meinhard #	OEM P/No.	Description
TRP-50-A0.5		Meinhard Low-uptake Type-A+ Concentric Glass Nebulizer, high sensitivity
TQP-50-A0.5		Meinhard Low-uptake Type-A+ Quartz Concentric Nebulizer, high sensitivity
SilQ+		High-Purity Quartz concentric nebulizer, with PFA fittings as standard, for ultra-trace analyses
TR-50-C1	47929001	Type C nebulizer, glass, 1 mL/min and 1 L/min @ 50 psi
TR-30-C2	47929002	Type C nebulizer, glass, 2 mL/min and 1 L/min @ 30 psi
TR-30-C3	47929003	Type C nebulizer, glass, 3 mL/min and 1 L/min @ 30 psi
TR-50-C0.5	47929004	Type C nebulizer, glass, 0.5 mL/min and 1 L/min @ 50 psi
TR-30-K3	47929005	Type K nebulizer, glass, 3 mL/min and 0.7 L/min @ 30 psi
TR-50-C3	47929006	Type C nebulizer, glass, 3 mL/min and 1 L/min @ 50 psi
F0		PFA Capillary with Micro Tubing. L=70cm, ID=0.25mm. Small-bore fitting for low uptakes.
F1		Detachable PFA capillary and 3-mm ID connector for general application, L=70cm, ID=0.8mm.
F2		Standard PFA Integrated Capillary, for self-aspiration or pumped samples, L=70cm, ID=0,50mm.
Fit Kit 3T		PTFE sample fitting and capillary for 4-mm OD (large-bore) liquid connection tubes
FG		PFA Gas/Argon Quick Connect Fitting for Meinhard nebulizers. Zero contamination, easy to use.
PMM 4000		PEEK MiraMist Parallel-Path Nebulizer
PRL-01		'PEARL' Cleaning Device for all types of concentric glass nebulizers
PRG-01		'PERGO' argon humidifier, for all nebulizers, NO pressure vessel - uses membrane permeation
ML116054	31026996	JY Spray Chamber with ground-glass Joint (also 31031326)
ML116055	31030809	Cyclonic Spray Chamber (also 31031279), top socket, 4 mm drain
2015-015NB	31-031-279	JY Cylonic Spray Chamber, no baffle
MM151030	31-023-549	PTFE Endcap for CG Nebulizer, 35/6 mm adaptor
ML116054	31026996	Spray Chamber Scott Type (also 31031326)
ML116055	31030809	Cyclonic Spray Chamber (also 31031279), top socket, 4 mm drain
ML116055ES		ML116055 with Easy-Seal threaded nebulizer inlet and Teflon PFA cap
ML116056		Cyclonic Spray Chamber, 50 mL, top socket, ground glass tapered (14/20) drain
ML116057	47929046	ML116055 with water jacket for cooling
ML116060	47929049	Cyclonic Spray Chamber, 20 mL, top socket, 4mm drain, Unifit inlet
ML116062	31031479	Cyclonic Spray Chamber, Axial, Taper to Top Socket, O-ring Inlet, 4 mm Drain
ES-1952-9001	ES-1952-9001	JY Demountable Torch by ESI, with Sapphire Injector ID = 1,5mm
MP16040	21-356-110	JY Fully Demountable Torch Assembly
ML116050	31023722	Quartz Torch Tube, Outer Dem 71mm
ML116051		JY FDT Quartz Outer Tube L=82.5mm
ML116052	31023723	Quartz Torch Tube, Inner, Demountable
ML116053	31023724	Teflon Centering Ring
MA1589	31021589	Alumina Injector, 2.0mm
ML116059		Quartz Injector, ID=2.0mm
ML116063	21356040	Pyrex Sheath Gas Tube
ML116070	31031099	Quartz Torch Tube, Outer, Panorama, L=105mm
ML135005	49010100	Quartz Torch, Fixed Injector, 1.5mm
MM116200	46702525	O-ring, Viton, Pack of 4
MM116210		O-ring Kit, Demountable torch (9 o-rings)
Tube-PVCSD-30	3710027200	PVC 2-Stop Tubing 0.030" Black/Black
Tube-PVCSD-40	3710034400	PVC 2-Stop Tubing 0.040" White/White
Tube-PVCSD-51	3710034500	PVC 2-Stop Tubing 0.051" Gray/Gray
Tube-PVCSD-65	3170034600	PVC 2-Stop Tubing 0.065" Blue/Blue
Tube-SolFSD-30	3710034800	SolventFlex 2-Stop Tubing 0.030" Black/Black
Tube-SolFSD-40	3710035000	SolventFlex 2-stop Tubing 0.040" (1.02 mm), White/White
Tube-SolFSD-51	3710035200	SolventFlex 2-Stop Tubing 0.051" Gray/Gray
Tube-SolFSD-65	3710067900	SolventFlex 2-stop Tubing 0.065" (1.65 mm), Blue/Blue

--	--	--