

<b>Parts list for Teledyne Leeman ICP instruments</b>		
<b>Meinhard P/N°</b>	<b>OEM P/No.</b>	<b>Descriptions</b>
SilQ+		High-Purity Quartz concentric nebulizer, with PFA fittings as standard, for ultra-trace analyses
TQ-30-C1		Concentric quartz nebulizer, Type C, with PFA argon and sample line connectors, high solids
TR-30-C2		Type C nebulizer, glass, 2 mL/min, high solids, replacement for # 120-00474-1
TR-30-K3		Type K nebulizer, glass, 3 mL/min and 1 L/min @ 30 psi, 4mm OD nebulizer sample tube
TR-30-K2		Type K nebulizer, glass, 2 mL/min and 1 L/min @ 30 psi, 4mm OD nebulizer sample tube
TQ-30-A2		Quartz Type-A concentric nebulizer, 2 mL/min, 1 L/min, 30 psi, high sensitivity
F0		PFA Capillary with Micro Tubing. L=70cm, ID=0.25mm. Small-bore fitting for low uptakes.
F1		Detachable PFA capillary and 3-mm ID connector for general application, L=70cm, ID=0.8mm.
F2		Standard PFA Integrated Capillary, for self-aspiration or pumped samples, L=70cm, ID=0,50mm.
Fit Kit 3T		PTFE sample fitting and capillary for 4-mm OD (large-bore) liquid connection tubes
FG		PFA Gas/Argon Quick Connect Fitting for Meinhard nebulizers. Zero contamination, easy to use.
PRL-01		'PEARL' Cleaning Device for concentric glass nebulizers - suits all types concentric glass nebs
PRG-01		'PERGO' argon humidifier - no pressure vessel - suits all nebulizer types
ML155025	120-3751	Spray Chamber, Radial, Scott Type, requires nebulizer adapter
ML155053	120-00475-2	Spray Chamber, Radial, Cyclonic, O-rings for 6mm nebulizer
ML155054	120-00475-1	Spray Chamber, Radial, Cyclonic, Baffle, O-rings for 6mm nebulizer
ML155055	318-00049	Spray Chamber, Axial, Scott Type (requires adapter cap; also 120-00201-1)
ML155058	120-00461-2	Spray Chamber, Axial, Cyclonic, for Hildebrand Grid
ML155059	120-00461-1	Spray Chamber, Axial, Cyclonic, with Baffle, for Hildebrand Grid (also 318-00092)
MM155100		Tubing (PVC) and quick disconnects for Leeman torch
ML155062	318-00062	Quartz Torch, one-piece Axial, AQ (187mm length), 2.5mm injector
ML155064	318-00064	Quartz Torch, one-piece Axial, AQ-aqueous, (160mm length)
ML155065	318-00065	Quartz Torch, Dual View with slot, AQ-aqueous
MP155066	318-00066	Purge tube, Axial, for Axial and DV Profile
ML155066	120-00496-1	Quartz Torch, Dual View, DURA
MP155067	318-00067	Purge tube, Radial, for DV Profile
MP155068	318-00068	Purge tube, Radial for Profile
ML155069	318-00069	Quartz Torch, Demountable - Axial, Direct-Reading Echelle
ML155070	318-00070	Quartz Torch, Demountable, Dual-View Direct-Reading Echelle
ML155020	120-00336-1	Quartz Torch, 1.5 mm injector (does not include tubing and fittings; 120-3748), radial view
ML155022	120-00336-2	Quartz Torch, 1.0 mm injector (does not include tubing and fittings; 120-3749), radial view
ML155023	120-00336-3	Quartz Torch, 2.5 mm injector, axial DRE (does not include tubing and fittings)
ML155043	318-0001	Quartz Torch tube, Original Demountable, radial view
Tube-PVCSD-20	309-00069-1	PVC tubing for the internal standard, orange/yellow tabs, ID=0.51mm
Tube-SolFSD-20	309-00069-2	SolventFlex tubing for the internal standard, orange/yellow tabs, ID=0.51mm
Tube-PVCSD-30	309-00069-4	PVC tubing for the drain, red/red tabs, ID=1.14mm
Tube-SolFSD-30	309-00069-5	SolventFlex tubing for the drain, red/red tabs, ID=1.14mm
Tube-PVCSD-45	309-00069-7	PVC tubing for the sample, black/black tabs, ID=0.76mm
Tube-SolFSD-45	309-00069-8	SolventFlex tubing for the sample, black/black tabs, ID=0.76mm