

Parts for Nu-instruments MC-ICP-MS		
Meinhard N°	OEM P/No.	Description
TRP-30-A0.5		Meinhard Type A+ Concentric Glass Nebulizer, with PFA argon and sample connectors
TQP-30-A0.5		Meinhard Type A+ Concentric Quartz Nebulizer, with PFA argon and sample connectors
SiIQ+		High-Purity Quartz concentric nebulizer, with PFA fittings as standard
TR-30-A1		Type A Concentric Glass Nebulizer, 1 mL/min, 1 L/min, 30 psi
TQ-30-A1		Quartz Type-A concentric nebulizer, 1 mL/min, 1 L/min, 30 psi
F0		PFA Capillary with encapsulated Micro Tubing. L=70cm, ID=0.25mm. For low uptakes.
F1		Detachable PFA capillary, for general application, L=70cm, ID=0.8mm.
F2		Standard PFA Integrated Capillary, for self-aspiration or pumped samples, L=70cm, ID=0,50mm.
FG		PFA Gas Quick Connect Fitting for Meinhard nebulizers. Zero contamination, easy to use.
PRL-01	8100-2051	'PEARL' Cleaning Device for concentric glass nebulizers - suits all types, all manufacturers
PRG-01		'PERGO' argon humidifier - keeps nebulizers running smoothly
PRG-01-MAPlus		'PERGO' argon high solids accessory including Meinhard Type A Plus high performance nebulizer
ES-4157-1000	ES-4157-1000	PC3 ESI <i>Peltier-cooled Quartz Spray Chamber</i> , with 12-Volt DC Power Supply
ML151016		Pyrex Mini-Cyclonic Spray Chamber, 20mL, Top Socket, 4mm drain, O-ring Inlet
ML151016Q		Quartz Mini-Cyclonic Spray Chamber, 20mL, Top Socket, 4mm drain, O-ring Inlet
ML170005	319-308	Std Quartz 1-Piece Torch with 1.5mm-ID Injector, O-ring ball joint
ML170025	319-311	Thermo PlasmaQuad / Nu Quartz Bonnet for 1-piece torch
ML170005R		Rebuild of ML170005 quartz torch
ML145005	3201192	Thermo PlasmaQuad/NU Quartz Torch, Axial, w/Ball Joint, 1.5mm Injector
ES-1223-0103-01		Complete Demountable ESI Quartz Torch, with 1-mm ID <i>quartz</i> injector in PFA base
ES-1223-0153-02		Complete Demountable ESI Quartz Torch, with 1,5-mm ID <i>quartz</i> injector in PFA base
ES-1223-0203-03		Complete Demountable ESI Quartz Torch, with 2-mm ID <i>quartz</i> injector in PFA base
ES-1203-0103-04		Complete Demountable ESI Quartz Torch, with 1-mm ID <i>sapphire</i> injector in PFA base
ES-1203-0183-06		Complete Demountable ESI Quartz Torch, with 1,8-mm ID <i>sapphire</i> injector in PFA base
ES-1203-0000		O-ring free Sapphire Injector, with PFA base, ID = 1,0 or 1,8 mm
ES-1202-1203		Complete torch, O-ring free Sapphire Injector with PFA base, ID = 1,5mm, RS12/5 adaptor
ES-1213-0183-07		Complete Demountable ESI Quartz Torch, with 1,8-mm ID <i>platinum</i> injector in PFA base
ES-1213-0203-08		Complete Demountable ESI Quartz Torch, with 2-mm ID <i>platinum</i> injector in PFA base
ES-1213-0000		O-ring free Platinum Injector, with PFA base, ID = 1,8 or 2,0 mm
ES-5509		PFA adaptor, from RS13/5 socket to ESI O-ring free torch base
ES-5510		PFA adaptor, from ESI spray chamber to RS13/5 torch base
ES-3000-8201	319-285	Nickel Sampler Cone
ES-3000-8202	319-497	Nickel Skimmer Cone
ES-3000-8206		Platinum Sampler Cone
ES-3000-8208		Platinum Skimmer Cone