

EPOND Price List, PE Optima 3x00 DV / XL and RL		
Meinhard P/No.	OEM P/No.	Description
TRP-50-A0.5		Meinhard Low-uptake Type-A+ Concentric Glass Nebulizer, high sensitivity
TQP-50-A0.5		Meinhard Low-uptake Type-A+ Quartz Concentric Nebulizer, high sensitivity
SilIQ+		High-Purity Quartz concentric nebulizer, with PFA fittings as standard, for ultra-trace analyses
TQ-30-A3	WE024371	Type A nebulizer, Quartz, 3 mL/min and 1 L/min @ 30 psi
TR-30-A3	00472020	Type A nebulizer, glass, 3 mL/min and 1 L/min @ 30 psi
TR-50-C1	0047-2022	Type C nebulizer, glass, 1 mL/min and 1 L/min @ 50 psi
TR-30-K3	N068-1574	Type K nebulizer, glass, 3 mL/min and 0.7 L/min @ 30 psi
TR-30-C3	N810-2011	Type C nebulizer, glass, 3 mL/min and 1 L/min @ 30 psi
TR-30-C0.5	N810-2012	Type C nebulizer, glass, 0.5 mL/min and 1 L/min @ 30 psi
Fit Kit 3	WE024372	Solution Inlet Adapter for 4 mm nebulizers, includes 1.2 m ETFE tubing.
F0		PFA Capillary with encapsulated Micro Tubing. L=70cm, ID=0.25mm. For low uptakes.
F1		Detachable PFA capillary, for general application, L=70cm, ID=0.8mm.
F2		Standard PFA Integrated Capillary, for self-aspiration or pumped samples, L=70cm, ID=0.50mm.
FG		PFA Gas Quick Connect Fitting for Meinhard nebulizers. Zero contamination, easy to use.
ES-4040-28	N0777219	PolyPro-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube
ES-2040-28	N0777049	PFA-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube
PMM 4000	N0775330	Burgener PEEK MiraMist Parallel-Path Nebulizer, 0,2 to 2,5 mL/min.
Snapper 1.251	N0777460	Nebulizer gas tubing clamp
MM21045		O-ring kit for GemTip Nebulizer, pack of 4
PRL-01	8100-2051	'PEARL' Cleaning Device for concentric glass nebulizers - suits all types, all manufacturers
PRG-01		'PERGO' argon humidifier - keeps nebulizers running smoothly
PRG-01-MAPlus		'PERGO' argon high solids accessory including Meinhard Type A Plus high performance nebulizer
ML130019XL	N812-2188	Cyclonic Glass Spray Chamber, Optima 3000DV, 3000XL
ML130024	218080256	Nebulizer Adapter - unifit for 6 mm o.d. nebulizer to 16 mm Unifit inlet
ML130019		Cyclonic Glass Spray Chamber, Optima 3000 RL, SCR
ES-4157-1000		PC3 Peltier-thermostatted Quartz Spray Chamber, with HF-res Nebulizer
ES-2044-0002		Threaded Drain line for spray chamber, PharMed, 0,75mm x 1m.
ML131001		Quartz Torch Tube, 3 slots, 138mm, Optima 3000 RL
ML131004	N0690568	Quartz Torch Tube, 1 slot, 138mm, Optima 3000 RL
ML131005		Quartz Torch Tube, no slot, 125mm, Optima 3000 RL
ML131012		Optima 3x00DV Quartz Torch, 2-slot
ML131014	NO69-1662	Optima 3x00DV Quartz Torch Tube, 1-slot
ML131015		Optima 3x00DV Quartz Torch Tube, 3-slot
ML131031	N0691664	Quartz Bonnet for 3000 Axial/DV
ML131032	N0695456	Quartz Bonnet for 3000 (after 9/21/94)
ML131033	N0695295	Quartz Bonnet for 3000 RL (before 9/21/94)
MP690780	N0690780	Optima 3000DV Injector Cyclonic Spray Chamber Support Adaptor
MM690769	N0690769	O-Ring Kit for Optima 3000XL/DV, packet of 9
ML132040	N0695441	Quartz Injector, 2.0mm, 3x00 XL and 3000 RL
ML132041	N0695442	Quartz Injector, 2.0 mm, 3x00 XL / DV, L=83mm
ML132044	N0695449	Quartz Injector, 1.6mm, 3x00 XL / DV, L=83mm
ML130046	N0695460	Quartz Injector, 1.2 mm, Optima 3000 XL/DV, 4300V, 5300V
MA5362	N069-5362	Alumina Injector, Optima 3x00 XL/DV ID=2.0, L=83mm
MA2186	N0582186	Alumina Injector, Optima 3x00 RL, ID=2.0. Also N8126039
MP2731	9992731	Axial purge window, Optima 2x00/3x00XL/3000SCX
MP6549	B0066549	Axial purge window, Optima 3000 DV, 4000, 5000, 7000 DV
MP1497	N0581497	Quartz Radial purge window, Optima 3000R, 3000 SCR, 3000 RL, 4x00V, 5x00V
MP0672	N0690672	Quartz Radial purge window, Optima 3000 DV - short tube
MP1689	N0691689	Quartz Axial purge window, Optima 3000 DV - long tube
Tube-PVCSD-30	09908587	Two-stop peristaltic pump tubing in PVC, black/black, pack of 12
Tube-PVCSD-30-F		Flared two-stop peristaltic pump tubing in PVC, black/black, 12,
Tube-PVCSD-45	N0680375	Two-stop peristaltic pump tubing in PVC, red/red, pack of 12
Tube-SolFSD-30	00473550	Two-stop peristaltic pump tubing in SolventFlex, black/black, pack of 12
Tube-SolFSD-30-F		Flared two-stop peristaltic pump tubing in SolventFlex, black/black, 12
Tube-SolFSD-45	N0680376	Two-stop peristaltic pump tubing in SolventFlex, red/red, pack of 12