

Shimadzu ICP-OES		
Meinhard #	OEM P/No.	Component Descriptions
TRP-50-A0.5		Meinhard Low-uptake Type-A+ Concentric Glass Nebulizer, high sensitivity
TQP-50-A0.5		Meinhard Low-uptake Type-A+ Quartz Concentric Nebulizer, high sensitivity
SilQ+		High-Purity Quartz concentric nebulizer, with PFA fittings as standard, for ultra-trace analyses
TR-30-C1	046-00092-02	'Meinhard' Concentric glass nebulizer, Type C, high solids
TQ-30-C1		'Meinhard' Concentric Quartz nebulizer, Type C, high solids
TR-30-A1		'Meinhard' Concentric glass nebulizer, Type A, high sensitivity
ES-4040-28		PolyPro-ST HF-resistant concentric Nebulizer, with detachable uptake tube
ES-2040-28		PFA-ST HF-resistant concentric Nebulizer, with detachable uptake tube
F0		PFA Capillary with Micro Tubing. L=70cm, ID=0.25mm. Small-bore fitting for low uptakes.
F1		Detachable PFA capillary and 3-mm ID connector for general application, L=70cm, ID=0.8mm.
F2		Standard PFA Integrated Capillary, for self-aspiration or pumped samples, L=70cm, ID=0,50mm.
Fit Kit 3T		PTFE sample fitting and capillary for 4-mm OD (large-bore) liquid connection tubes
FG		PFA Gas/Argon Quick Connect Fitting for Meinhard nebulizers. Zero contamination, easy to use.
PMM-4000		PEEK MiraMist Parallel-Path Nebulizer, for high solids and slurry samples
PRL-01		'PEARL' Cleaning Device for concentric glass nebulizers – suits all types concentric glass nebs
PRG-01		'PERGO' argon humidifier - no pressure vessel - suits all nebulizer types
	NOTE :	Meinhard spray chambers require pumped drain
ML151026		Pyrex Cyclonic Spray Chamber, 50 mL, Top Socket, Inner Baffle, 4 mm Drain, O-rings
ML151026ES		ML151026 with Easy-Seal threaded nebulizer inlet and Teflon PFA cap
ML151026Q		Quartz Cyclonic Spray Chamber, 50 mL, Top Socket, Inner Baffle, 4 mm Drain, O-rings
ML151027		Pyrex Cyclonic Spray Chamber, 50 mL, Top Socket, No Baffle, 4 mm Drian, O-rings
ML151027ES		ML151027 with Easy-Seal threaded nebulizer inlet and Teflon PFA cap
ML151027Q		Quartz Cyclonic Spray Chamber, 50 mL, Top Socket, No Baffle, 4 mm Drain, O-rings
ME41571000	ES-4157-1000	PC3 Peltier-cooled Cyclonic Spray Chamber, with 12-Volt Power Supply
ME91005700	ES-9100-5700	PC3-LT Peltier-cooled Quartz Spray Chamber, reaches -20 degrees C.
2016-0005	20470272	Quartz Torch, radial, with 1,5mm injector ID, #ML160005
ML160008	21181448	Quartz Mini-Torch for ICPE-9000, L = 117mm., # 211-81448SSG
ML160007	21146605	Shimadzu Quartz Torch, Axial, with curved Inlet to Tapered Injector
Tube-PVCSD-25		PVC 2-stop tubing, 0.025"/0.64mm, orange/white
Tube-PVCSD-25-F		PVC 2-stop tubing, 0.25"/0.64mm, orange/white, flared both ends
Tube-PVCSD-40		PVC 2-stop tubing, 0.040"/1.02mm, white/white
Tube-PVCSD-35-3		PVC 3-stop tubing, 0.035"/0.89mm, orange/orange
Tube-PVCSD-35-3F		PVC 3-stop tubing, 0.035"/0.89mm, orange/orange, flared both ends
Tube-PVCSD-45-3		PVC 3-stop tubing, 0.045"/1.14mm, red/red