

<b>Spectro SOP ICP-OES instruments</b>		
<b>Meinhard #</b>	<b>OEM P/No.</b>	<b>Component Descriptions</b>
SilQ+		High-Purity Quartz concentric nebulizer, PFA fittings as standard, for ultra-trace analyses
TQ-30-C1		Concentric quartz nebulizer, Type C, with PFA argon and sample line connectors
TR-30-A1		Type A concentric glass nebulizer, 1 mL/min and 1 L/min @ 30 psi
TR-30-K2		Type K concentric glass nebulizer, 2 mL/min and 1 L/min @ 30 psi
TR-30-C1	76060510	Type C high-solids concentric glass nebulizer, 1 mL/min and 1 L/min @ 30 psi
PMM 4000		PEEK MiraMist nebulizer
MM175100		Liquid inlet fitting and argon tubing with quick disconnect for TR-30-C1
ES-4040-28		PolyPro-ST HF-resistant concentric nebulizer, with detachable sample capillary
ES-2040-28		PFA-ST HF-resistant concentric nebulizer, with detachable sample capillary
ES-2030-28		PFA-HS High-Solids HF-res concentric nebulizer, with detachable sample capillary
8100-2051		'PEARL' Cleaning Device for concentric glass nebulizers - for all manufacturer's types
PERGO		'PERGO' argon humidifier - no pressure vessel - suits all nebulizer types
ML175010		Spray Chamber, Slurry
ML175012	48105003	Spray Chamber, Upward Conical Open
ML175014	48105006	Spray Chamber, Large Double Pass
ML175018	48106008	Spray Chamber, Large, Double Pass, Jacketed, Threaded Ports
ML175026	48105001	Cyclonic Spray Chamber, Top Socket, Pumped Drain
ML175026ES		ML175026 with Easy-Seal threaded nebulizer inlet and Teflon PFA Cap
ML175026Q		Quartz Cyclonic Spray Chamber, Socket
ML175026QJ	USA00850	Quartz Cyclonic Spray Chamber, Jacketed, Threaded Ports
ES-2256-0470		47-mm diameter PFA spray chamber for ICP-OES, with 'PureCap' end-cap
ES-3167-1001		Quartz cyclonic spray chamber, can be used with PC3 Peltier cooling unit
ES-4157-1000		PC3 Peltier-cooled Quartz Spray Chamber, with 12-Volt DC Power Supply
ML175005	75060596	Spectro Quartz Torch, 1.8 Injector, suits Spectroflame
ML175006	75060545	Quartz Torch, Radial, Demountable, 80mm Length
ML175007	48205001	Quartz torch, 1.0mm Injector, for organics
ML175008		Quartz Torch, Radial, with 2.5mm injector
ML175009	75360521	Quartz torch, organic (does not include gas fittings)
ML175043	48206002	Quartz Torch, SOP, Demountable, No Flare, 100mm Length
MM175055		Tubing and gas connectors 44302072, 44302077
ML175021	48105008	Quartz Injector, 1.5 mm
ML175022	48105010	Quartz Injector, 1.8mm
ML175023	48105009	Quartz Injector, 2.5mm
Tube-PVCSD-15-F		PVC 2-Stop Tubing 0.015"/0.38mm Orange/Green, Flared both ends, pkt 12
Tube-PVCSD-35-F		PVC 2-Stop Tubing 0.035"/0.89mm Orange/Orange, Flared both ends, pkt 12
Tube-PVCSD-51-F		PVC 2-Stop Tubing 0.051"/1.3mm Gray/Gray, Flared both ends
Tube-PVCSD-51		PVC 2-Stop Tubing 0.051"/1.3mm Gray/Gray, pkt 12
Tube-PVCSD-69		PVC 2-Stop Tubing 0.069"/1.75mm Blue/Green, pkt 12