

<b>Meinhard accessories for Thermo ICP-OES instruments</b>		
<b>Meinhard #</b>	<b>OEM P/No.</b>	<b>ICAP-6300, 6500, DUO</b>
TRP-50-A0.5		Meinhard Low-uptake Type-A+ Concentric Glass Nebulizer, high sensitivity
TQP-50-A0.5		Meinhard Low-uptake Type-A+ Quartz Concentric Nebulizer, high sensitivity
SilQ+		High-Purity Quartz concentric nebulizer, with PFA fittings as standard, for ultra-trace analyses
TR-30-K2	842312051431	Type K nebulizer, glass, 2 mL/min and 1 L/min @ 30 psi
TR-50-A2		Type A nebulizer, glass, 2 mL/min and 1 L/min @ 50 psi
TR-50-C2	842312051331	Type C nebulizer, glass, 2 mL/min and 1 L/min @ 50 psi, high solids
MP080963	824312051321	V-Groove Nebulizer, Noordermeer, 6 mm O.D., 30 psi
F0		PFA Capillary with Micro Tubing. L=70cm, ID=0.25mm. Small-bore fitting for low uptakes.
F1		Detachable PFA capillary and 3-mm ID connector for general application, L=70cm, ID=0.8mm.
F2		Standard PFA Integrated Capillary, for self-aspiration or pumped samples, L=70cm, ID=0,50mm.
Fit Kit 3T		PTFE sample fitting and capillary for 4-mm OD (large-bore) liquid connection tubes
FG		PFA Gas/Argon Quick Connect Fitting for Meinhard nebulizers. Zero contamination, easy to use.
ES-4040-28		PolyPro-ST HF-resistant concentric Nebulizer, with detachable uptake tube
ES-2040-28		PFA-ST HF-resistant concentric Nebulizer, with detachable uptake tube
PMM 4000	842312051471	MiraMist Parallel-Path Nebulizer in PEEK. 1 L/min @ 45 psi
PRL-01	8100-2051	'PEARL' Cleaning Device for concentric glass nebulizers - suits all types, all manufacturers
PRG-01		'PERGO' argon humidifier - keeps nebulizers running smoothly
PRG-01-MAPLus		'PERGO' argon high solids accessory including Meinhard Type A Plus high performance nebulizer
ML128010	842312051411	Glass Cyclonic Spray Chamber, iCAP 6x00
ML128010ES		ML128010 with Easy-Seal threaded nebulizer inlet and Teflon PFA cap
ML128012	842312051311	Glass Cyclonic Spray Chamber, with Baffle, iCAP 6x00
ML128012ES		ML128012 with Easy-Seal threaded nebulizer inlet and Teflon PFA cap
ML128021	842312051251	Spray Chamber right-angled Adapter, for iCAP 6x00 Duo
ML128020	842312051151	Spray Chamber straight Adapter, for iCAP 6x00 Duo
ML128007	842312051241	Duo Quartz Torch, iCAP 6x00, older type
ML128040	842312051341	Quartz Center Tube, 1.0 mm, iCAP 6x00
ML128041	842312051351	Quartz Center Tube, 1.5 mm, iCAP 6x00
ML128042	842312051371	Quartz Center Tube, 2.0 mm, iCAP 6x00, L=97mm
ML128050	842312051841	EMT Duo Quartz Torch, iCAP 6x00
ML128070	842312051941	Quartz Center Tube, 1.0 mm, iCAP 6x00
ML128071	842312051951	Quartz Center Tube, 1.5 mm, iCAP 6x00
ML128072	842312051971	Quartz Center Tube, 2.0 mm, for iCAP 6x00 EMT torch
ML128073	842312051591	Radial View Window/Sleeve, for iCAP Duo,
Tube-PVCSD-25	842312051511	PVC 2-stop tubing, 0.025"/0.64mm, orange/white
Tube-PVCSD-25-F		PVC 2-stop tubing, 0.25"/0.64mm, orange/white, flared both ends
Tube-PVCSD-40	842312051521	PVC 2-stop tubing, 0.040"/1.02mm, white/white
Tube-PVCSD-35-3	13017700	PVC 3-stop tubing, 0.035"/0.89mm, orange/orange
Tube-PVCSD-35-3F		PVC 3-stop tubing, 0.035"/0.89mm, orange/orange, flared both ends
Tube-PVCSD-45-3	13107701	PVC 3-stop tubing, 0.045"/1.14mm, red/red
MP500QD		Internal liquid fitting for spray chamber drain