

Thermo X-series ICP-MS		
Meinhard #	OEM P/No.	Part Descriptions
TRP-50-A0.5		Meinhard Low-uptake Type-A+ Concentric Glass Nebulizer, replacement for G3266-65003
TQP-50-A0.5		Meinhard Low-uptake Type-A+ Quartz Concentric Nebulizer, upgrade for G3266-65003
SiIQ+		High-Purity Quartz concentric nebulizer, with PFA fittings as standard
TR-30-A3	1040894	Type A concentric glass nebulizer, 3 mL/min and 1 L/min @ 30 psi
TQ-30-A3		Type A concentric Quartz nebulizer, 3 mL/min and 1 L/min @ 30 psi
TR-30-C3		Type C nebulizer, glass, 3 mL/min and 1 L/min @ 30 psi
TR-30-A1	1201318	Type A concentric glass nebulizer, 1 mL/min and 1 L/min @ 30 psi
TQ-50-A0.5		Type A concentric quartz nebulizer, 0.5mL/min, low uptake
F0		PFA Capillary with encapsulated Micro Tubing. L=70cm, ID=0.25mm. For low uptakes.
F1		Detachable PFA capillary, for general application, L=70cm, ID=0.8mm.
F2		Standard PFA Integrated Capillary, for self-aspiration or pumped samples, L=70cm, ID=0,50mm.
FG		PFA Gas Quick Connect Fitting for Meinhard nebulizers. Zero contamination, easy to use.
PMM 4000		Burgener PEEK MiraMist Parallel-Path Nebulizer
PRL-01		'PEARL' Cleaning Device for concentric glass nebulizers
PRG-01		Pergo argon high solids accessory without nebulizer. Includes gas line connector kit.
PRG-01-MAPlus		Pergo argon high solids accessory including Meinhard Type A Plus high performance nebulizer
ML145016	3200841	Quartz Spray Chamber, Scott type, Jacketed for coolant
ML145032	3200727	Quartz Coupling, Ball/Socket, 90 deg, for 3200841 Spray Chamber
ML145022	3204313	Glass Spray Chamber, Conical, Jacketed
ML145023	3600374	Quartz Spray Chamber, Conical, Jacketed
ML145024	3600170	Quartz Spray Chamber, Conical, Peltier Cooled, X-Series
ML145025	3600169	Quartz Coupling, Ball/Socket, 17 deg
ML145026	4600508	Quartz Conical Spray Chamber, Peltier-cooled, Small Drain, X-Series 2
ML145027	3600093	Cyclonic Spray Chamber, Quartz, X-Series
MP500QD		Internal liquid fitting for spray chamber drain
ES-4190-60-29		PC3 Sample Introduction Kit for XSeries2 with Quartz Spray Chamber and PFA-ST nebulizer
ES-9190-60-29		PC3-LT Sample Introduction Kit for XSeries2 with Quartz Spray Chamber and PFA-ST nebulizer
ML145050	3601145	Standard Quartz Torch for X-series ICP-MS, 2 locating beads, 1.5mm Injector, ES-1203-0015
ML145051	3601146	Demountable Quartz Torch, ICPMS, X-Series 2, with pin
ME12030020	3601147	Standard Quartz Torch for X-series ICP-MS, 2 locating beads, 2mm Injector
ME12021203	ES-1202-1203	HF-res Semi-Demountable Torch with Sapphire Injector, O-ring free
ML145053	3601149	Quartz Bonnet, X series 2 , ES-1202-9001
ML145010	3600970	Quartz Torch, ICPMS X Series, Demountable, 14 mm
ML145011	3600969	Quartz Torch, ICPMS X Series, w/Ball Joint, 1.5mm Injector
ML145012	3600791	Quartz Torch, ICPMS X Series, w/Ball Joint, 2.0mm Injector
ME12020301	ES-1202-0301	O-ring free, Precision Quartz Torch, X-series, ESI
ES-3000-1235	3601219	PlasmaScreen Silver Screen for torch
ES-4397-1038		PVC 2-Stop Tubing, flared ends, ID=0.38 mm, Orange/Green
ES-4397-1051		PVC 2-Stop Tubing, flared ends, 0.51 mm, Orange/Yellow
Tube-PVCSD-51		PVC 2-Stop Tubing 0.051" Gray/Gray
Tube-PVCSD-65		PVC 2-Stop Tubing 0.065" Blue/Blue
Tube-PVC-10-FA		PVC 2-stop (72 mm) Tubing, 0.010" (0.254 mm), Orange/Blue, Flared Ends, 19 cm
Tube-PVC-075-FA		PVC 2-stop (72 mm) Tubing, 0.0075" (0.19 mm), Orange/Red, Flared Ends, 19 cm