

VG PQ series, Axiom, ExCell, PlasmaTrace		
EPOND N°	OEM P/No.	Part Descriptions
TRP-50-A0.5		Meinhard Type A+ Concentric Glass Nebulizer, with PFA argon and sample connectors
TQP-50-A0.5		Meinhard Type A+ Concentric Quartz Nebulizer, with PFA argon and sample connectors
SilQ+		High-Purity Quartz concentric nebulizer, with PFA fittings as standard
TR-30-K2		K-Type Concentric Glass Nebulizer, 2 ml/min.
TR-30-K1	1040894	K-High Flow, AX Conc. Glass Nebulizer, 1 ml/min.
TR-30-C1		C-Type standard Concentric Glass Nebulizer, 1 ml/min.
PMM 4000		PEEK MiraMist Parallel-Path Nebulizer
PRL-01	8100-2051	'PEARL' Cleaning Device for concentric glass nebulizers - suits all types, all manufacturers
PRG-01		'PERGO' argon humidifier - keeps nebulizers running smoothly
PRG-01-MAPlus		'PERGO' argon high solids accessory including Meinhard Type A Plus high performance nebulizer
ML145022	3204313	VG Jacketed Pyrex Spray Chamber with Impact Sphere, for PQ/X-series
ML145023	3600374	VG Quartz Jacketed Spray Chamber, with Impact Sphere, for PQ/X-series
ML145024	3600170	VG Quartz Conical PQ/X-series Spray Chamber, for Peltier system
ML145016	3200841	Quartz Coaxial VG Spray Chamber, Scott type, Jacketed for Coolant
EP-2124-0351	ES-2124-0351	35mm PFA spray chamber and PureCap End Cap for PQ2 and PQ3
EP-6124-0001	ES-6124-0001	Spray Chamber Cooling Block, anodized Al, plus Mounting Bracket No. 6124-1002
EP-6124-1002	ES-6124-1002	Mounting Bracket for attachment of Al cooling block to PQ2
EP-6124-2002	ES-6124-2002	Mounting Bracket for attachment of Al cooling block to VG Axiom
ML145005	3201192	VG Standard Quartz Torch, 1.5mm Injector, with Ball-Joint
ML170005	319-308	Quartz Torch for PQ and Nu-instruments, 1.5mm Injector, Ball Joint with O-ring seal
ML145006	3201732	VG Standard Quartz Torch with Ball Joint, 2mm Injector
ML145008	1200080	VG ICPMS Demountable Quartz Torch Body
ME12020201	ES-1202-0201	O-ring Free Precision Quartz Torch, ESI
ML170025	319-311	Quartz Bonnet for PQ and Nu-instruments
ML145029	3203789	VG Torch Shield, for PlasmaScreen option, suits VG PQ ExCell, also P/N. 3600884
ML145025	3600169	VG Quartz Elbow
ML145032	3200727	VG Elbow, Ball/Socket Quartz
ES-4398-1076		Peri-pump tubes for Sample, 2-tags, black-black, ID = 0,76 mm, pkt 12
ES-4398-1102		Peri-pump tubes for Sample, 2-tags, white-white, ID = 1,02 mm, pkt 12
ES-4398-1165		Peri-pump tubes for Drain, 2-tags, blue-blue, ID = 1,65 mm, pkt 12
ES-4397-1038		PVC 2-Stop Tubing, flared ends, ID=0.38 mm, Orange/Green, pkt 12
ES-4397-1051		PVC 2-Stop Tubing, flared ends, 0.51 mm, Orange/Yellow, pkt 12
Tube-PVCSD-51		PVC 2-Stop Tubing 0.051" Gray/Gray, pkt 12
Tube-PVCSD-65		PVC 2-Stop Tubing 0.065" Blue/Blue, pkt 12