

Agilent / Varian Axial Vista, 720-ES and axial-view Liberty		
Meinhard #	OEM P/No.	Descriptions
TRP-50-A0.5		Meinhard Low-uptake Type-A+ Concentric Glass Nebulizer, high sensitivity
TQP-50-A0.5		Meinhard Low-uptake Type-A+ Quartz Concentric Nebulizer, high sensitivity
SiIQ+		High-Purity Quartz concentric nebulizer, with PFA fittings as standard, for ultra-trace analyses
TQ-30-C1		Concentric quartz nebulizer, Type C, with PFA argon and sample line connectors, high solids
TR-30-C2	2010076500	Type C nebulizer, glass, 2 mL/min and 1 L/min @ 30 psi, high solids, 4mm OD sample tube
TR-30-K3	2010081600	Type K nebulizer, glass, 3 mL/min and 1 L/min @ 30 psi, 4mm OD nebulizer sample tube
TR-30-K2	2010069500	Type K nebulizer, glass, 2 mL/min and 1 L/min @ 30 psi, 4mm OD nebulizer sample tube
TQ-30-A2		Quartz Type-A concentric nebulizer, 2 mL/min, 1 L/min, 30 psi, high sensitivity
F0		PFA Capillary with Micro Tubing. L=70cm, ID=0.25mm. Small-bore fitting for low uptakes.
F1		Detachable PFA capillary and 3-mm ID connector for general application, L=70cm, ID=0.8mm.
F2		Standard PFA Integrated Capillary, for self-aspiration or pumped samples, L=70cm, ID=0.50mm.
Fit Kit 3T		PTFE sample fitting and capillary for 4-mm OD (large-bore) liquid connection tubes
FG		PFA Gas/Argon Quick Connect Fitting for Meinhard nebulizers. Zero contamination, easy to use.
PMM-4000	G3161-80000	PEEK MiraMist Parallel-Path Nebulizer, for high solids and slurry samples
8100-2051		'PEARL' Cleaning Device for concentric glass nebulizers – suits all types concentric glass nebs
PERGO		'PERGO' argon humidifier - no pressure vessel - suits all nebulizer types
ML180019	7910043700	Pyrex Cyclonic Spray Chamber, Baffle, 8 mm Top, Unifit
ML180019ES		ML180019 with Easy-Seal threaded nebulizer inlet and Teflon PFA cap
ML180021	2010081700	Pyrex Cyclonic Spray Chamber, 8 mm Top, O-ring Inlet
ML180021ES		ML180021 with Easy-Seal threaded nebulizer inlet and Teflon PFA cap
ML180022		Cyclonic glass spray chamber, with Baffle, O-ring seal for nebulizer
ML180013		Pyrex Cyclonic Spray Chamber, with Baffle, Top Socket, Unifit
ML130024	218080256	Nebulizer Adapter - Unifit for 6 mm o.d. nebulizer to 16 mm inlet
ML180023	9910070400	Cyclonic glass spray chamber, with Coolant Jacket
ME21280471	ES-2128-0471	47 mm PFA Spray Chamber and ESI 'PureCap' end cap
ME21281000	ES-2128-1000	Mounting Bracket for PFA Spray Chamber
ML180017	2010090400	Quartz Torch, ICP AX, w/Right Angle End, 2.3 mm injector (also 2010092800)
ML180050	2010094800	Quartz Torch, ICP AX, High Solids, with 90-degree bend
ML180014	2010081200	Quartz Torch, ICP AX, 2.3 mm injector
ML180047	2010077500	Quartz Injector, 0.8mm ID, for organic samples
ML180015	2010083400	Quartz Torch, Outer Tubes for Demountable torch
MA6759	2410067590	Alumina Injector 1.8 mmID x 4mmOD x 129mm
ML180025	2010070700	Quartz Bonnet, ICP (also 2010070790)
ML180010	2010072900	Quartz Torch, Outer Tubes for Demountable torch
ML180046	9910057000	Quartz Injector, 1.4mm ID
MP180048	9910057300	Injector Support Adapter PTFE
ME18430180	ES-1843-0180	ESI Demountable Torch sapphire injector, ID = 1,8mm, O-ring free
ME18423020	ES-1842-3020	ESI Demountable Torch platinum injector, ID = 2mm, O-ring free
ME1842	ES-1842	ESI Demountable Torch Quartz Tubes, for use with O-ring free injector base
ME43992300	ES-4399-2300	Male Barbed PFA connector for Transfer Line
ME4399200251	ES-4399-2002-51	PFA Elbow
ME18421100	ES-1842-1100	Complete Demountable Torch, with O-ring free injector base and Sapphire Injector
Tube-PVCSD-30	3710027200	PVC 2-Stop Tubing 0.030" Black/Black
Tube-PVCSD-40	3710034400	PVC 2-Stop Tubing 0.040" White/White
Tube-PVCSD-51	3710034500	PVC 2-Stop Tubing 0.051" Gray/Gray
Tube-PVCSD-65	3170034600	PVC 2-Stop Tubing 0.065" Blue/Blue
Tube- SolFSD-30	3710034800	SolventFlex 2-Stop Tubing 0.030" Black/Black
Tube-SolFSD-40	3710035000	SolventFlex 2-stop Tubing 0.040" (1.02 mm), White/White
Tube-SolFSD-51	3710035200	SolventFlex 2-Stop Tubing 0.051" Gray/Gray
Tube-SolFSD-65	3710067900	SolventFlex 2-stop Tubing 0.065" (1.65 mm), Blue/Blue