

Agilent / Varian radial Vista, 725-ES and radial-view Liberty		
Meinhard #	OEM P/No.	Descriptions
TRP-50-A0.5		Meinhard Low-uptake Type-A+ Concentric Glass Nebulizer, high sensitivity
TQP-50-A0.5		Meinhard Low-uptake Type-A+ Quartz Concentric Nebulizer, high sensitivity
SIQ+		High-Purity Quartz concentric nebulizer, with PFA fittings as standard, for ultra-trace analyses
TQ-30-C1		Concentric quartz nebulizer, Type C, with PFA argon and sample line connectors, high solids
TR-30-C2	2010076500	Type C nebulizer, glass, 2 mL/min and 1 L/min @ 30 psi, high solids, 4mm OD sample tube
TR-30-K3	2010081600	Type K nebulizer, glass, 3 mL/min and 1 L/min @ 30 psi, 4mm OD nebulizer sample tube
TR-30-K2	2010069500	Type K nebulizer, glass, 2 mL/min and 1 L/min @ 30 psi, 4mm OD nebulizer sample tube
TQ-30-A2		Quartz Type-A concentric nebulizer, 2 mL/min, 1 L/min, 30 psi, high sensitivity
F0		PFA Capillary with Micro Tubing. L=70cm, ID=0.25mm. Small-bore fitting for low uptakes.
F1		Detachable PFA capillary and 3-mm ID connector for general application, L=70cm, ID=0.8mm.
F2		Standard PFA Integrated Capillary, for self-aspiration or pumped samples, L=70cm, ID=0.50mm.
Fit Kit 3T		PTFE sample fitting and capillary for 4-mm OD (large-bore) liquid connection tubes
FG		PFA Gas/Argon Quick Connect Fitting for Meinhard nebulizers. Zero contamination, easy to use.
PMM-4000	G3161-80000	PEEK MiraMist Parallel-Path Nebulizer, for high solids and slurry samples
PRL-01		'PEARL' Cleaning Device for concentric glass nebulizers – suits all types concentric glass nebs
PRG-01		'PERGO' argon humidifier - no pressure vessel - suits all nebulizer types
ML180019	7910043700	Pyrex Cyclonic Spray Chamber, Baffle, 8 mm Top, Unifit
ML180019ES		ML180019 with Easy-Seal threaded nebulizer inlet and Teflon PFA cap
ML180021	2010081700	Pyrex Cyclonic Spray Chamber, 8 mm Top, O-ring Inlet
ML180021ES		ML180021 with Easy-Seal threaded nebulizer inlet and Teflon PFA cap
ML180022		Cyclonic glass spray chamber, with Baffle, O-ring seal for nebulizer
ML180013		Pyrex Cyclonic Spray Chamber, with Baffle, Top Socket, Unifit
ML130024	218080256	Nebulizer Adapter - UniFit for 6 mm o.d. nebulizer to 16 mm inlet
ML180023	9910070400	Cyclonic glass spray chamber, with Coolant Jacket
ME21280471	ES-2128-0471	47 mm PFA Spray Chamber and ESI 'PureCap' end cap
ME21281000	ES-2128-1000	Mounting Bracket for PFA Spray Chamber
ML180005	2010069690	Varian radial-view ICPQuartz Torch, Low Flow, 1.5mm Injector, also # 2010069600
ML180005B	2010069610	Varian radial-view ICP Quartz Torch, Low Flow, 1.5mm Injector, with ball joint
ML180028	2010117400	Varian radial-view ICP Torch, Injector 0.8mm
ML180029	2010117500	Quartz Torch, ICP Radial, 1.4 mm injector
ML180025	2010070790	Varian Radial Quartz Torch Bonnet, also # 20-100707-00, ML180025
ML180046	2010073000	Varian Radial Liberty Quartz Injector 1.8mm ID
MA6759	2410067500	Varian Radial Liberty Alumina Injector 1.8mm ID, also # 241007590
MM180048	9910057300	Varian Radial Liberty Teflon Injector Holder for demountable torch
Tube-PVCSD-30	3710027200	PVC 2-Stop Tubing 0.030" Black/Black
Tube-PVCSD-40	3710034400	PVC 2-Stop Tubing 0.040" White/White
Tube-PVCSD-51	3710034500	PVC 2-Stop Tubing 0.051" Gray/Gray
Tube-PVCSD-65	3170034600	PVC 2-Stop Tubing 0.065" Blue/Blue
Tube-SolFSD-30	3710034800	SolventFlex 2-Stop Tubing 0.030" Black/Black
Tube-SolFSD-40	3710035000	SolventFlex 2-stop Tubing 0.040" (1.02 mm), White/White
Tube-SolFSD-51	3710035200	SolventFlex 2-Stop Tubing 0.051" Gray/Gray
Tube-SolFSD-65	3710067900	SolventFlex 2-stop Tubing 0.065" (1.65 mm), Blue/Blue